



Ankle-Arm Blood Pressure

Participant Id#:

Acrostic:

**AABPDT3**  
 Date:   /   /      
Month Day Year

**Pulse Obliteration Pressure:**

*Inflate the cuff slowly until the pulse is no longer audible. Inflate to 30 mm Hg above the level at which pulse sound disappeared.*

**1 Systolic Readings:**

Record in this order

	<b>Systolic</b>	
<b>RBRACH3</b> Right brachial		<input type="text"/> <input type="text"/> <input type="text"/>
<b>RDPEIS3</b> Right dorsalis pedis		<input type="text"/> <input type="text"/> <input type="text"/>
<b>RPTIB3</b> Right posterior tibial		<input type="text"/> <input type="text"/> <input type="text"/>
	<b>Systolic</b>	
<b>LDPEDI</b> Left dorsalis pedis		<input type="text"/> <input type="text"/> <input type="text"/>
<b>LPTIB3</b> Left posterior tibial		<input type="text"/> <input type="text"/> <input type="text"/>
<b>LBRACH3</b> Left brachial		<input type="text"/> <input type="text"/> <input type="text"/>

**2 Procedure was: AABPCMP3**

- Completed successfully → Skip to Technician ID
- Not completed → Continue to #3

**3 Reason procedure was not completed:**

<b>RIGHT SIDE:</b>	<b>YES</b>	<b>NO</b>
Unable to occlude	<input type="radio"/>	<input type="radio"/> <b>RLUOCCL3</b>
Ulceration	<input type="radio"/>	<input type="radio"/> <b>RLULCER3</b>
Amputation	<input type="radio"/>	<input type="radio"/> <b>RLAMPUT3</b>
Unable to locate DP	<input type="radio"/>	<input type="radio"/> <b>RULOCDP3</b>
Unable to locate PT	<input type="radio"/>	<input type="radio"/> <b>RULOCPT3</b>

Other:

<b>LEFTSIDE:</b>	<b>YES</b>	<b>NO</b>
Unable to occlude	<input type="radio"/>	<input type="radio"/> <b>LLUOCCL3</b>
Ulceration	<input type="radio"/>	<input type="radio"/> <b>LLULCER3</b>
Amputation	<input type="radio"/>	<input type="radio"/> <b>LLAMPUT3</b>
Unable to locate DP	<input type="radio"/>	<input type="radio"/> <b>LULOCDP3</b>
Unable to locate PT	<input type="radio"/>	<input type="radio"/> <b>LULOCPT3</b>

Other:

**AABPTID3** Technician ID#:

**AABPRID3** Reviewer ID#:

**AABPDID3** Data Entry ID #: